

## **Homedics Tanda Luxe Anti Ageing Skin Rejuvenation Photofacial device**

### **When can I expect to see results?**

---

When treating the signs of aging with the Luxe, immediate and long term results can be seen.

During treatment:

Users will feel a warm sensation along with a gentle vibration on the skin that some users have describe as calming. Many users experience an increase in hydration and skin plumping immediately.

After 1 week of treatment:

Most users will start to see an improvement in skin texture while the majority of users will see an increase in the healthy feel of their skin and an increase in skin smoothness.

After 2 weeks of treatment:

Many users see a reduction in the overall signs of aging, particularly an increase in skin firmness.

After 4 weeks of treatment:

Many users see an overall better skin condition compared to before treatment. This includes an improvement in the healthy appearance of the skin and a reduction in fine lines and wrinkles.

### **What is light therapy?**

---

Light therapy is the application of light energy for therapeutic benefits. It is a proven technology that has been used in clinics for rehabilitation and aesthetic applications for over 35 years.

### **Are there any harmful side effects?**

---

There are no known harmful side effects of light therapy. Thousands of scientific studies have proven this technology to be a safe and highly effective way to treat a variety of conditions. It is recommended to do a light sensitive test as outlined in the user manual before use to make sure you are not sensitive to light.

### **Can I use my normal creams and lotions?**

---

Yes, apply your creams and lotions after treatment with the Luxe device as your skin is now better prepared to absorb them

---

### **How long does the Luxe treatment head last for?**

The Luxe treatment heads last for 300 hours which is the equivalent of 6,000 three-minute cycles.

---

### **Can it treat through makeup or clothing?**

No. Laser and light devices cannot penetrate clothing sufficiently to provide an effective treatment. For best results cleanse your skin before treatment.

---

### **Will the LED's used in the Luxe burn my skin?**

No. Under normal operating circumstances, the temperature remains at a comfortable level. If the temperature exceeds this level, a temperature sensor in the device terminates the treatment cycle.

---

### **Is there any UV danger when using the Luxe?**

No. The Luxe does not use UV light. Luxe uses 660nm red LED light.

---

### **Do I have to wear the goggles?**

The goggles have been provided for your convenience and comfort as some people find the light bright. Company sponsored testing showed that, when used as directed, the device did not pose risk of damage to the eyes.

---

### **How long will a full charge last?**

A fully charged battery will provide approximately 8-10, three-minute cycles, a total of 24-30 minutes.

---

### **Can I take my Luxe when I travel?**

---

### **Is this device only for individual use?**

The benefit of the Luxe is that it can be shared by the whole family. Individual treatment heads may be purchased separately. The treatment heads can be easily cleaned between uses with a damp cloth and mild detergent.

---

### **Are the treatment heads interchangeable?**

---

Yes, the hand piece for both the Clear+ and Luxe are the same. Therefore you are able to interchange treatment heads on the same system offering several skin care solutions in one.

### **How long will the device itself last?**

---

The treatment heads last for 6,000 treatments; however the hand piece itself depends on usage and care.

### **What all comes included in my Luxe Kit?**

---

Your Luxe Kit will come with the following items:

- Luxe cordless rechargeable handpiece
- Luxe treatment head
- Recharging stand
- Power supply
- Goggles
- User manual

### **Can I use the device if I am pregnant or nursing?**

---

No, do not use the Luxe if you are pregnant or nursing.

### **Is there a maximum amount of treatment time per day?**

---

It is recommended to use the Luxe once a day for the best results, however there is no maximum on the amount of treatment time per day and no harm in over treating.

### **Is it alright if I overlap treatment areas?**

---

Yes, it is absolutely fine if your treatment areas overlap.

### **What is the difference between the Prestige and Professional models?**

---

Our Tanda Prestige system is our smaller model that offers 20 LED's in the treatment head, while our Professional model is about 1 and a half times larger and offers 36 LED's in the treatment head. The main difference is the surface area of the treatment head, with the Professional model it allows for a reduction in treatment time.

### **Is the 30 day money back guarantee applicable to all Tanda products?**

---

The 30 day money back guarantee is applicable on all products sold through [www.tanda.com](http://www.tanda.com) and is only valid if purchased direct from this site. If purchased from another reseller, their return policy would apply.

### **Can I use both the Luxe and Clear+ treatments in the same day? Is there any down time required between treatments?**

---

Yes, you can absolutely use both the Luxe and Clear+ in the same day. In fact, since there is no down time needed, you can do them back to back.

### **What is the difference between the Luxe and the older model – Tanda Regenerate**

---

Luxe is the next generation of the breakthrough Tanda Regenerate anti-aging device . While both use 660nm red LED to increase the production of collagen and elastin, the new Luxe is now 5 times more powerful and offers 2 new technologies – Sonic Vibration and gentle warming to maximize results.