

Homedics Tanda Clear Professional Acne clearing system

When can I expect to see results?

You can see a difference in your complexion in just 24 hours!

Does the Clear+ Prevent breakouts from occurring?

Yes! The Clear+ uses blue LED light that is scientifically proven to kill P. acnes bacteria, treat existing blemishes and prevent further breakouts. If used as part of your daily skincare routine, you can prevent your breakouts from ever happening.

What should I expect to feel or see after using the Clear+?

First 24-48 Hours:

During treatment users will feel a warm sensation on the skin that some have described as calming. Many users experience a decrease in blemish size and an improvement in their complexion.

After 1 week of treatment:

Most users will see a decrease in the number of blemishes. The majority will also see a reduction in redness and skin oiliness, and will continue to see an overall improvement in skin condition. When used regularly, Clear+ significantly reduces the amount of acne-causing bacteria on the skin helping to prevent future breakouts.

Is the Clear+ a one-spot treatment?

No, the surface area of the Clear+ is 2.59in² /16.7cm², so it can be used to treat several acne affected areas at once. The larger treatment surface also helps to prevent further breakouts in the surrounding area. For a spot treatment see theTanda ZAP.

What is light therapy?

Light therapy is the application of light energy for therapeutic benefits. It is a proven technology that has been used in clinics for rehabilitation and aesthetic applications for over 35 years.

Are there any harmful side effects?

There are no known harmful side effects of light therapy. Thousands of scientific studies have proven this technology to be a safe and highly effective way to treat a variety of conditions.

Can I use my normal creams and lotions?

Yes, your creams and lotions can be applied before or after treatment.

How long does the Clear+ treatment head last for?

The Clear+ treatment heads last for 300 hours which is the equivalent of 6,000 three-minute cycles.

Can it treat through makeup or clothing?

No. Laser and light devices cannot penetrate clothing sufficiently to provide an effective treatment.

Will the LED's used in the Clear+ burn my skin?

No. Under normal operating circumstances, the temperature remains at a comfortable level. If the temperature exceeds this level, a temperature sensor in the device terminates the treatment cycle.

Is there any UV danger when using the Clear+?

No. The Clear+ does not use UV light. The Clear+ uses 414nm blue LED light.

Do I have to wear the goggles?

The goggles have been provided for your convenience and comfort as some people find the light bright. Company sponsored testing showed that, when used as directed, the device did not pose risk of damage to the eyes.

How long will a full charge last?

A fully charged battery will provide approximately 8-10, three-minute cycles, a total of 24-30 minutes.

Can I take my Clear+ when I travel?

Yes, the Clear+'s light weight design and compact form make traveling easy. The Clear+ also features a dual 110-220V transformer and a plug converter so you can use the device all over the world.

Is this device only for individual use?

The benefit of the Clear+ is that it can be shared by the whole family. Individual treatment heads may be purchased separately. The treatment heads can be easily cleaned between uses with a damp cloth and mild detergent.

Are the treatment heads interchangeable?

Yes, the hand piece for both the Clear+ and Luxe are the same. Therefore you are able to interchange treatment heads on the same system offering several skin care solutions in one. The new clear+ and Luxe treatment heads are also compatible with our old prestige devices.

How long will the device itself last?

The treatment heads last for 6,000 treatments; however the hand piece itself depends on usage and care.

What all comes included in my clear+ Kit?

Your Tanda Clear+ Kit will come with the following items:

- Clear+ cordless rechargeable handpiece
- Clear+ treatment head
- Recharging stand
- Power supply
- Goggles
- User manual

Can I use the device if I am pregnant or nursing?

No, do not use the Clear+ if you are pregnant or nursing.

I have deeper, cystic acne, is the Clear+ right for me?

The Clear+ is recommended for mild to moderate acne. It is not recommend for deeper, more cystic acne. Please see the chart below to find out if the Clear+ is right for you:

Acne Type	Description of Condition	Is Clear+ Right for you?
Mild Acne	<ul style="list-style-type: none">• Several inflamed pimples – red in color• Less than 20 whiteheads, or less than 15 inflammatory red pimples or less than 30 total pimples not all inflamed	Yes
Moderate Acne	<ul style="list-style-type: none">• Any inflamed pimples (red in color) and pustules (visible accumulation of pus in skin)• 20 to 100 whiteheads, or 15 to 50 inflammatory pimples, or 30 to 125 total pimples not all inflamed	Yes

Severe Nodular Acne	<ul style="list-style-type: none"> Inflamed pimples and pustules (visible accumulation of pus in skin) with a few deep nodular lesions (solid mass can be felt under skin – can sometimes be raised). Greater than 5 cysts (solid mass of skin like a knot, can be raised or felt under the skin), or total whiteheads count Greater than 100, or total inflammatory count Greater than 50, or Greater than 125 total lesions 	No, consult your dermatologist for the best prescribed treatment.
Severe Cystic Acne	<ul style="list-style-type: none"> Many nodular cystic lesions (with signs of scarring) 	No, consult your dermatologist for the best prescribed treatment.

Is there a maximum amount of treatment time per day?

It is recommended to use the Clear+ twice a day for the best results, however there is no maximum on the amount of treatment time per day and no harm in over treating.

Is it alright if I overlap treatment areas?

Yes, it is absolutely fine if your treatment areas overlap.

Is the 30 day money back guarantee applicable to all Tanda products?

The 30 day money back guarantee is applicable on all products sold through www.tanda.com and is only valid if purchased direct from this site. If purchased from another reseller, their return policy would apply.

Can I use both the Clear+ and Luxe treatments in the same day? Is there any down time required between treatments?

Yes, you can absolutely use both the Clear+ and Luxe in the same day. In fact, since there is no down time needed, you can do them back to back.

What is the difference between the Clear+ and the older model - Tanda Clear?

The difference between the two are the technologies. While both use 414nm blue LED to kill the bacteria causing acne, the new Clear+ introduces 2 new technologies – Sonic Vibration and Heat to open up the pores so the bacteria have no where to hide. The new Clear+ head is compatible with the older devices.