



Condom facts

Lubricated Odour masker Easy-on shape Extra lubrication Security Sensitivity Pleasure Latex free Benefits

		Lubricated	Odour masker	Easy-on shape	Extra lubrication	Security	Sensitivity	Pleasure	Latex free	Benefits
Security	Extra Safe									Slightly thicker and extra lubricated to provide extra reassurance
	Comfort XL									Extra long and wider, providing better comfort for those who need it
	Close Fit									A tighter fit for extra confidence
Sensitivity	Elite									Super thin and extra lubricated for greater sensitivity
	Fetherlite									Ultra fine for enhanced sensitivity
Pleasure	Tingle									Creates a gentle tingling sensation, heightening stimulation for both partners
	Pleasuremax									Uniquely placed ribs and dots plus shaped and flared end for ultimate pleasure for both partners
	Sensation									500 raised dots for enhanced pleasure and intensity
	Select Flavours									Coloured and flavoured for extra excitement
	Ribbed									Ribbed for greater stimulation
Special	Performa									Contains a special benzocaine cream inside the teat end to help control climax and prolong sexual excitement, made to make you last longer
	Avanti Ultima									Made from polyisoprene, Avanti Ultima offers people with latex allergies a softer and more versatile option, compared with traditional materials




Size and shape facts

Shaped for extra pleasure

Special fit

Latex free Easy-on



All measurements are nominal and are lay flat width.
 a = 30mm from bead b = 70mm from bead
 c = 120mm from bead d = widest part of the bulbous section
 Thickness measured in microns.

The illustrations shown are for size and shape comparison use, to demonstrate the difference between the varying types of condoms. The dimensions indicated are estimated for guidance only.

Extra Safe
Thickness 75

Fetherlite
Thickness 55

Performa
Thickness 65

Select Flavours
Thickness 70

Tingle
Thickness 65

For further information about Durex condoms and lubricants:

- Telephone: **0800 338 739**
- E-mail: **info@durex.co.uk**
- Visit: **www.durex.co.uk**
- Write to: **FREEPOST DUREX**

