

Mouth Spray

PAIN RELIEF OF MOUTH ULCERS AND SMALL ORAL LESIONS

INDICATION

Aloclair Plus® Mouth Spray has a mechanical action (barrier effect) indicated for the management of the symptoms caused by aphthous ulcers, oral stomatitis, mouth sores, oral irritations and small oral injuries due to braces and dentures. It is particularly indicated for mouth ulcers and lesions difficult to reach.

COMPOSITION

AQUA, GLYCERIN, HYDROGENATED STARCH HYDROLYSATE, SODIUM CITRATE, POLYSORBATE 80, PVM/MA COPOLYMER, MALTODEXTRIN, PENTYLENE GLYCOL, GLUCONOLACTONE, PEG-40 HYDROGENATED CASTOR OIL, PVP, SODIUM HYDROXIDE, SODIUM BENZOATE, ALOE BARBADENSIS LEAF JUICE, GLYCERYL CAPRYLATE/CAPRATE, POLYCARBOPHIL, XANTHAN GUM, ALTHAEA OFFICINALIS ROOT EXTRACT, SODIUM HYALURONATE, SUCRALOSE, AROMA, CALCIUM GLUCONATE, POTASSIUM SORBATE.

HOW TO USE

For single spread lesions. When using the product for the first time take the top off and apply vaporizer. After correct oral hygiene, direct vaporizer towards the lesion inside the mouth and spray 2-4 times to cover the affected area and to create barrier film. Do not touch the lesion with the tongue for at least 2 minutes to allow the formation of the protective film. The vaporizer allows the product to cover the oral lesions with maximum efficiency. Apply the product 3-4 times a day or as needed, until the symptoms disappear. After thirty days of use, suspend the application for at least one week, before restarting the treatment. For best result, avoid eating and drinking for at least one hour after application.

WARNINGS

Carefully close the cap after each use. There are no known side effects deriving from the use of the product, however it is contra-indicated in patients with a known history of hypersensitivity to any of the ingredients. If the symptoms persist, or in case of symptoms worsening, ask for a medical advice.

The product is not intended to be ingested, but topically applied in the oral cavity. If the product is swallowed accidentally, no adverse effects are anticipated. Do not use after the expiry date, which is indicated in the pack. Keep out of the sight and reach of children. Do not use if the package is damaged or open. The product should be used in a single patient to avoid cross contamination. Avoid contact with eyes.

The product has not been evaluated during pregnancy and breastfeeding. If you are pregnant/breastfeeding, consult your doctor before use. Suitable for children from 4 years of age.

STORAGE

The expiry date refers to the intact product correctly stored. Store in the original pack, in a dry place away from light and heat source. Do not refrigerate the product. Store at room temperature.

Size: 15ml \bigcirc bottle with enveloped vaporizer.

Biokosmes s.r.l. – socio unico

Via dei Livelli 1, 23842 Bosisio Parini (LC) ITALY

Distributed by: Alliance Pharmaceuticals Ltd

Avonbridge House, Bath Road, Chippenham, Wiltshire SN15 2BB, United Kingdom

Aloclair®, Alliance and associated devices are trademarks of the Alliance Pharma PLC group.

Symbol	1	C € 0477	Ţ	\subseteq	LOT
Meaning	legal manufacturer of medical device	CE mark from notified body	Consult the instruction for use for important information	Expiry date	Batch number



Rev 09-2022 TR2646127 0ALL2FI00145